



Recombinant Putative carbonic anhydrase 5 (cah-5)

Product Code	CSB-MP608802CXY
Abbreviation	cah-5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10462
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	GPGSDHGWGY DENNGPDTWQ GKCQNHLKQS PIDIRAPDVD YALLHRMHFL NYDMDGKIEL SNTGRTL FAG GFESWQHKQP MIQGGGLKHR YKLAQFHLHW GQNDAVGSEH AMGSLHYP AE LHLVHVREGL TLKEALSRPD GLAVVGVFLA KTNDPVANKF SPISERLHDL RHSGNKTELK NFRTKYVLPL DTEAFYRYEG SLTTPDCSEA VIWTVLAEPM AISSHQLHLL RQLHNKELVK SDKNYRPLQP LNGRRIQYRP SKLDRAMICS SFATTSVIIT YVAILSAMLI
Source	Mammalian cell
Target Names	cah-5
Protein Names	Recommended name: Putative carbonic anhydrase 5 EC= 4.2.1.1 Alternative name(s): Carbonate dehydratase 5
Expression Region	21-310
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.