



Recombinant RNA-binding protein rnp-1 (rnp-1)

Product Code	CSB-EP608824CXY-B
Abbreviation	rnp-1
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10667
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MPSKLVGNL PDNVDSNKLK QVFQPFCKVT ECDIVKNYAF VHIEEDDVPD IITRLTGYTI DGKVVNIKKS TSKLRPTPGM PNRCFRCQSD EHRTQPQPD PTNNQKTENG VQTLKFDLTS GAGVKRSAGD PIIDSAKRIA YGAQSVVEPE IPQPMDDLQ ALYQEYQLSR QRYVYYRDRL LKEMEAKQHG STAGFALSSS STVPVPVASA PPGATQLSAA PVSYPNAPP VIASINAPYA VASNLRAPYA LQSAPYASAA SAPYGSVTPA GAPSNVMTTQ QYLQQIQHQQ ATGSPAPVPA PPRLY
Source	E.coli
Target Names	rnp-1
Protein Names	Recommended name: RNA-binding protein rnp-1
Expression Region	1-305
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.