



Recombinant *Oryza sativa* subsp. *japonica* Germin-like protein 3-7 (GER7)

Product Code	CSB-BP608864OFG
Abbreviation	GER7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10BU2
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	DPEPI QDFCVAVPRA GGEASPAYPG FPCKPASAVV SDDFFFAGLA AAGSTDNPFQ ASLKPGNVEA FPALNTLGVA INRVDLAPGG VNPLHSHPRA AELVHVITGR MLVGFVSTAG KYYSKVVGEG ETFAIPRGLM HFQYNPGNAS ARAMTVFNSQ LPGVVPAATA LFGADPEIPD AVLAKSFQVD AEIIKLLKSK FKK
Source	Baculovirus
Target Names	GER7
Protein Names	Recommended name: Germin-like protein 3-7 Alternative name(s): Germin-like protein 7 Short name= OsGER7
Expression Region	36-233
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.