



Recombinant *Oryza sativa* subsp. *japonica* Putative B3 domain-containing protein Os03g0619850 (Os03g0619850, LOC_Os03g42250)

Product Code	CSB-MP608870OFG
Abbreviation	Os03g0619850, LOC_Os03g42250
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10GP0
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAGGGSRMKK SCACCKRYLE HLGKMSCFL IRMTTDSMHS MIIPDRFVNH FGGKIPGTIK LESPNGILYV VEVTECMNKT VLQCGWEAFV DAHHIKVGDS LLFRHIENSC FEVMILDSDG SENVSLKSNR NGVSDSQES EDSEGPAGPP YILSWKSKSR LSSLQKKIHK EKVRSIQSEV PIYVAIMNKS NIGLTSSPCQ LELGARYAAA VHLPDRRQAV VLQRGAAMGH RDADQERQMH HQAVPDQWLE QDS
Source	Mammalian cell
Target Names	Os03g0619850
Protein Names	Recommended name: Putative B3 domain-containing protein Os03g0619850
Expression Region	1-253
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.