



Recombinant *Oryza sativa* subsp. japonica Protein IN2-1 homolog A (Os03g0283200, LOC_Os03g17480)

Product Code	CSB-EP608879OFG
Abbreviation	Os03g0283200, LOC_Os03g17480
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10N44
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAAAAAPRSS GKEALPAALG SASEPPRLFD GTTRLYICYF CPFAQRAWII RNFKGLQDKI ELVGIDLQDK PAWYKEKVEE QGTVPSELEHN GKIMGESLDL IKYIDSHFEG PALLPEDPEK RQFADELIAY ANAFTKALYS PLISKADLSA ETVAALDKIE AALSKFGDGP FFLGQFSLVD IAYVTIIERI QIYYSHIRKY EITNGRPNLE KFIEEINRIE AYTQTKNDPL YLLDLAKTHL KVA
Source	E.coli
Target Names	Os03g0283200
Protein Names	Recommended name: Protein IN2-1 homolog A
Expression Region	1-243
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.