



Recombinant *Oryza sativa* subsp. japonica Probable protein phosphatase 2C 28 (Os03g0137200, LOC_Os03g04430)

Product Code	CSB-BP608889OFG
Abbreviation	Os03g0137200, LOC_Os03g04430
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q10S32
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MLAAVMDYFR SCWGPRSPAG HRVRGSDVAG RQDGLLWYKD AGQLVTGEFS MAVVQANNLL EDQSQVESGA LSMAEPGPQG TFIGVYDGHG GPETARFIND HMFHHLRRFA TEHKCMSTDV IRKAFQATEE GFLSLVSKQW SLKPQIAAVG SCCLVGVICS GTLYVANLGD SRAVLGRFVK STGEVVATQL SSEHNACYEE VRQELQASHP DDPQIVVLKH NVWRVKGLIQ ISRSIGDVYL KRPEYNREPL HSKFRLRETF KRPILSSEPA IAVHQIQPND HFVIFASDGL WEHLSNQEAV DLVQNNPRNG IARRLVKVAM QEAACKREMR YSDLKKIDRG VRRHFHDDIT VIVVFLDSNA ISKANWSRGP SVSLRGGGVT LPANSLAPFS TPTVLSSTY
Source	Baculovirus
Target Names	Os03g0137200
Protein Names	Recommended name: Probable protein phosphatase 2C 28 Short name= OsPP2C28 EC= 3.1.3.16
Expression Region	1-399
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.