



Recombinant Uncharacterized calcium-binding protein B0563.7 (B0563.7)

Product Code	CSB-BP608924CXY
Abbreviation	B0563.7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q11083
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MRSAKVGVAR QLETKKPQTG RKISTSSRGT IHSQQSQPED IQMKYTRKEL KEYRQLFNMF DTDGSGAIGN EELKQAMISI GLHANKAEID NVIKEVDADG NGEIDFEEFC ACMKKSQIV KSTNEELIRE CFEIFDQDRN GIITENEFKY IAKEFGDFDD ELAEKVFREL DVSANGHLSA DQFATIVEDY LLNDPKHDID TGDS DVERYD DRHDDRASPM PNHLSTVPE
Source	Baculovirus
Target Names	B0563.7
Protein Names	Recommended name: Uncharacterized calcium-binding protein B0563.7
Expression Region	1-229
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.