



Recombinant Putative serine/threonine-protein kinase C01C4.3 (C01C4.3)

Product Code	CSB-YP608925CXY
Abbreviation	C01C4.3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q11090
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MTPEKTKEVF QMFKGAPTRR KSENAHSHRV VLNQLQNHPP RNATQSPQRQ PRTSESSMDF PRSALRRNST DTHVFANRTH NHRDINVPAC SSEEDRVSTA RRNSLFVKRG SVTMEPIKKV DLEEVYTVNK QLGTGRFGFI KLAEHKQSKQ RIAIKFFPRP QTKQADRVRE YNYSFFLSPH QNIIDTYEGM FQSSDDTAYF FVQEFPCPRAS LREAVEATNQ AGIGEANTKK VFAAVLSAIE FMHDENLVHR NLKAENILIF DANDYSKVKV TDFGLTRKVD TTVKYLEYVN NYHAAELCDT VVNEKLVVNK STDIWALGII FFYCMKGKFP WQKASIMCKP YWEWEQWLKR KNPALPKFN PFSEKALKLF KKSLTPRFKD RWTAKDMRKC LAKEKLLKSV KVAVPYY
Source	Yeast
Target Names	C01C4.3
Protein Names	Recommended name: Putative serine/threonine-protein kinase C01C4.3 EC=2.7.11.1
Expression Region	1-407
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.