



Recombinant Uncharacterized protein C03B1.2 (C03B1.2)

Product Code	CSB-BP608937CXY
Abbreviation	C03B1.2
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q11109
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	≥85% (SDS-PAGE)
Sequence	MRNLKNLSKM LLAGKISKYY RVKYLKFKCL VSVTTAEPNM EQLFQKYEEK LKLQNVQARI WEEKHTYIKS GKAVFKVNP PKQERTFRWV IEAPASSCQT SRGYSVKFTP SRVKRAKIP RMYNMKRCFQ LRRDSRIALS DAIFQPEKKT NTMDHLFTAS TAPEYSEKEL RQIEDFQSFF DNWDDDDYYTE QEKYPVIYDL ITIQISDCLD LMQQAVNGRI EVYQEVYKQT SNPIYLAKAV FIKFLDFDRL SPLHKTYLML SEECEFPNCY VEGILRRHYE
Source	Baculovirus
Target Names	C03B1.2
Protein Names	Recommended name: Uncharacterized protein C03B1.2
Expression Region	1-280
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.