



Recombinant Putative hydroxypyruvate isomerase (C05D11.5)

Product Code	CSB-EP608944CXY
Abbreviation	C05D11.5
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q11185
Product Type	Recombinant Protein
Immunogen Species	Caenorhabditis elegans
Purity	>85% (SDS-PAGE)
Sequence	MSGTNRVAAN LNMLFTNLPL LQRYGAAASA GFKLVEVSIP YTEPAEKLRE AADEYHLKHT LINAPPGNWD DGFRGLASLK SAKKEFRKSL DTAIEYAKAL GCCRVMHVMAG IPKSDDLEN AHQTYSENVF FAAEKFKCHK LICLIEPINK YTIPGYHLNN YEDAMDVIQM DQSNLKIYQY DTFHAQQING QIGAIMRKLK DYIGYIQVAQ VPNRGACDTR GEIDYHFIFD EIRSINPSWV IGAEYLDAKP SFWNIENMSL SF
Source	E.coli
Target Names	C05D11.5
Protein Names	Recommended name: Putative hydroxypyruvate isomerase EC= 5.3.1.22
Expression Region	1-262
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.