



Recombinant *Oryza sativa* subsp. japonica Double-stranded RNA-binding protein 3 (DRB3)

Product Code	CSB-YP609153OFG
Abbreviation	DRB3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0DJA3
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MKKKSAPTPL PPETANTSPA PIGATAGIRV ENCYVFKSRL QEYAQKAGLQ TPEYHTSKEG PSHEPVFKST VVINNTSYGS LPGFSNRKAA EQSAAEVALM EIVKSIPANA NIPAVQETGL CKNLLQEY AQ KMNYAIPSYI CTKPASGLAP FLCTVEIGGI QYIGAAARTK KDAEIKAAART ALLAIQQQSE GSANGATKYI VVP GKRVGKE VEKRPIETPK PLKAKKGGFK KKWNKRKFMK KDGQAVDVEK DEARVAGDAH DSDVLMQPTV ITQEASCRTL FLQPCEEAKR VEDEPPRDIE MVQPDKENQH SDAALVQPDD EARVEQEPSR DISVVQPNEE AISAKQEPSI DAATLQPKEE AMKTGCVALV LCLNISVDPP DVIKISPCVR KECWIDPFMS AAPKALETIG KTLHSQYERW QPKARYKLQL DPTLEEV
Source	Yeast
Target Names	DRB3
Protein Names	Recommended name: Double-stranded RNA-binding protein 3 Alternative name(s): dsRNA-binding protein 3
Expression Region	1-437
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.