



# Recombinant *Oryza sativa* subsp. *japonica* Probable protein phosphatase 2C 68 (Os09g0325700, LOC\_Os09g15670)

<b>Product Code</b>	CSB-EP609423OFG
<b>Abbreviation</b>	Os09g0325700, LOC_Os09g15670
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0J2L7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
<b>Purity</b>	≥85% (SDS-PAGE)
<b>Sequence</b>	MSMAEVCCDS AVVVGAEEA RARARAGRRR RAGVEGAGRW NATATAAGVA AEEAATRKRR ASGGEAGLVV VAKRHGAASV AGRRREMEDA VSLREAF AAP ANGEVAAARC DFYGVFDGHG CSHVADACRE RMHELVAEEM GAGSPAAAAR EPASWTETME RSFARMDAEV IAGCRAESGS CRCEGQKCDH VGSTAVVAVV EESRVVAVV GDSRAVLVCRG GAPVQLSSDH KPDRPDELER IEAAGGRVIF WEGARVLGVL AMSRSIGDAY LKPYPVTAVPE VTVTGRSDFD ECLILASDGL WDVVSNEAAC EVAQSCLRRG RQRWCAEAAA VLTKLALARR SSDNISVVVV DLRRGNAL
<b>Source</b>	E.coli
<b>Target Names</b>	Os09g0325700
<b>Protein Names</b>	Recommended name: Probable protein phosphatase 2C 68 Short name= OsPP2C68 EC= 3.1.3.16
<b>Expression Region</b>	1-358
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.