



Recombinant *Oryza sativa* subsp. *japonica* Probable protein phosphatase 2C 68 (Os09g0325700, LOC_Os09g15670)

Product Code	CSB-YP609423OFG
Abbreviation	Os09g0325700, LOC_Os09g15670
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0J2L7
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MSMAEVCCDS AVVVGAEAEA RARARAGRRR RAGVEGAGRW NATATAAGVA AEEAATRKRR ASGGEAGLVV VAKRHGAASV AGRRREMEDA VSLREAF AAP ANGEVAAARC DFIGVFDGHG CSHVADACRE RMHELVAEEM GAGSPAAAAR EPASWTETME RSFARMDAEV IAGCRAESGS CRCEGQKCDH VGSTAVVAVV EESRVVAVV GDSRAVLVCRG GAPVQLSSDH KPDRPDELER IEAAGGRVIF WEGARVLGVL AMSRSIGDAY LKPYYTAVPE VTVTGRSDFD ECLILASDGL WDVVSNEAAC EVAQSCLRRG RQRWCAEAAA VLTKLALARR SSDNISVVVV DLRRGNAL
Source	Yeast
Target Names	Os09g0325700
Protein Names	Recommended name: Probable protein phosphatase 2C 68 Short name= OsPP2C68 EC= 3.1.3.16
Expression Region	1-358
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.