



Recombinant *Oryza sativa* subsp. *japonica* Cinnamyl alcohol dehydrogenase 7 (CAD7)

Product Code	CSB-MP609427OFG
Abbreviation	CAD7
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0JA75
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MAPTTTATAA AEQAPPPQHT RKA VGLAAHD DSGHLTPIRI SRRKTGDDDDV AIKVLYCGIC HSDLHTIKNE WRNAVYPVVA GHEITGVVTE VGKNVARFKA GDEVGVGCMV NTCGGCESCR DGCENYCSGG VVFTYNSVDR DGTRTYGGYS DAVVVSQRFV VRFPS SAGGG AGAALPLDSG APLLCAGVTV YAPMRQHGLC EAGKHVGVVG LGGLGHVAVK FARAFGMRVT VISTSPVKRQ EALERLGADG FIVSTNASEM KAAMGTMHGI INTASASTSM HSYLALLKPK GKMILVGLPE KPLQIPTFAL VGGGKILAGS CMGSISETQE MIDFAAEHGV AADIELIGAD EVNTAMERLA KGDVRYRFVV DIGNTLRSD
Source	Mammalian cell
Target Names	CAD7
Protein Names	Recommended name: Cinnamyl alcohol dehydrogenase 7 Short name= OsCAD7 EC= 1.1.1.195 Alternative name(s): Protein FLEXIBLE CULM 1
Expression Region	1-379
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.