



Recombinant *Oryza sativa* subsp. japonica Homeobox-leucine zipper protein HOX17 (HOX17)

Product Code	CSB-MP609428OFG
Abbreviation	HOX17
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0JB92
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	≥85% (SDS-PAGE)
Sequence	MMERAEDLRL SLSLSSPLIA PRTHHVAMLF HAPPEKRFLE MPLLPAAKRS EVVAAEEERA GLRGGGGSDE EDGGCGIDGS RKKLRLSKDQ SAVLEDSFRE HPTLNPRQKA TLAQQLGLRP RQVEVWFQNR RARTKCLKQTE VDCEFLKRCC ETLTEENRRL QKEVQELRAL KLVSPHLYMN MSPPTTLTMC PSCERVSNTN NNSAAAAAD RRGIRTTTAA GGGSVVDTAA DGGILCHRPI AVRPPQS
Source	Mammalian cell
Target Names	HOX17
Protein Names	Recommended name: Homeobox-leucine zipper protein HOX17 Alternative name(s): HD-ZIP protein HOX17 Homeodomain transcription factor HOX17 OsHox17
Expression Region	1-247
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.