



Recombinant *Oryza sativa* subsp. *japonica* Cyclin-B1-3 (CYCB1-3)

Product Code	CSB-EP609437OFG-B
Abbreviation	CYCB1-3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0JNK6
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MASRRQWGAG AGAVDVPAQE RKGKKAATDH LLRAFFAVFF VLEGLSWRRR RRRDGEGVFL LRYVACLLCL QAAKIARRPK TTTVVAQQPP RIRRALADVS NLVNGRAALP VVNRQAAAA AADKCRKPIK QRNENNKAAK PEVIVISSDS EKHKKNPAQR AASRRAPIQT LTSILTKCSR ASDGVISPKK ELIIDIDASD SHNELAVVDY VEDIYRFYRN TENTYRPLCT YMVSQTEINE RMRAILDWL IEVHYRLMLM PETLYLTVYI IDQYLSLENV PRKELQLVGV SAMLIACKYE ETWAPLVKDF LVISDNSFSR QQVLSTEKSI LNKLQWNLTV PTMYMFILRY LKAALGDEEL EHMTFFYAEL ALVQYSMLFF APSVIAAAA YAARCTLGLS PLWSDLLEYH TGLAEPQLE CARRLVSLHA AAPESRQKVV YKKYASPKLG AVSLHSPAKK LLPPPSPVAA
Source	E.coli
Target Names	CYCB1-3
Protein Names	Recommended name: Cyclin-B1-3 Alternative name(s): CYCB1;1 G2/mitotic-specific cyclin-B1-3 Short name= CycB1;3
Expression Region	1-470
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.