



Recombinant *Oryza sativa* subsp. japonica Bowman-Birk type bran trypsin inhibitor (RBBI3.3)

Product Code	CSB-MP609440OFG
Abbreviation	RBBI3.3
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0JR25
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	ME RPWKCCDNK RLPTKPDPPQ WRCNDELEPS QCTAACKSCR EAPGPFPGKL ICEDIYWGAD GPGFCTPRPW GDCCDKAFCN KMNPTCRCM DEVKECADAC KDCQRVESSE PPRYVCKDRF TGHPGPVCKP R
Source	Mammalian cell
Target Names	RBBI3.3
Protein Names	Recommended name: Bowman-Birk type bran trypsin inhibitor Alternative name(s): OSE727A Protein RBBI3-3 RBTI
Expression Region	119-251
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.