



# Recombinant Gorilla gorilla gorilla NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 9 (NDUFB9)

|                          |   |
|--------------------------|---|
| <b>Product Code</b>      | CSB-MP609496GGZ   |
| <b>Abbreviation</b>      | NDUFB9  |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |
| <b>Uniprot No.</b>       | Q0MQE9  |
| <b>Product Type</b>      | Recombinant Protein   |
| <b>Immunogen Species</b> | Gorilla gorilla gorilla (Western lowland gorilla)   |
| <b>Purity</b>            | >85% (SDS-PAGE)   |
| <b>Sequence</b>          | AFLASGRYL THQQKVLRLY KRALRHLESW CVQRDKYRYF ACLMRARFEE HKNEKDMAKA TQLLKEAEEE FWYRQHPQPY IFPDSPPGGTS YERYDCYKVP EWCLDDWHPS EKAMYPDYFA KREQWKKLRR ESWEREVKQL QEETPPGGPL TEALPPARKE GDLPLWWYI VTRPRERPM   |
| <b>Source</b>            | Mammalian cell  |
| <b>Target Names</b>      | NDUFB9  |
| <b>Protein Names</b>     | Recommended name: NADH dehydrogenase [ubiquinone] 1 beta subcomplex subunit 9 Alternative name(s): Complex I-B22 Short name= CI-B22 NADH-ubiquinone oxidoreductase B22 subunit  |
| <b>Expression Region</b> | 2-179   |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.   |
| <b>Protein Length</b>    | Full Length of Mature Protein   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.   |