



Recombinant Escherichia coli O6:K15:H31 6-phosphofructokinase (pfkA)

Product Code	CSB-MP609768EGY
Abbreviation	pfkA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TAE8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MIKKIGVLTS GGDAPGMNAA IRGVRSALT EGLEVMGIYD GYLGLYEDRM VQLDRYSVSD MINRGGTFLG SARFPEFRDE NIRAVAIENL KKRIGDALVV IGGDGSYMGA MRLTEMGFPC IGLPGTIDND IKGTDYTIGF FTALSTVVEA IDRLRDTSSS HQRISVVEVM GRYCGDLTLA AAIAGGCEFV VVPEVEFSRE DLVNEIKAGI AKGKKHAIVA ITEHMCVDVE LAHFIEKETG RETRATVLGH IQRGGSPVPY DRILASRMGA YAILLLAGY GGRCVGIQNE QLVHHDIIIDA IENMKRPFKG DWLDCAKKLY
Source	Mammalian cell
Target Names	pfkA
Protein Names	Recommended name: 6-phosphofructokinase Short name= Phosphofructokinase EC= 2.7.1.11 Alternative name(s): Phosphohexokinase
Expression Region	1-320
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.