



Recombinant Escherichia coli O6:K15:H31 Rhamnulokinase (rhaB)

Product Code	CSB-BP609769EGY
Abbreviation	rhaB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TAF9
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MTFRNCVAVD LGASSGRVML ARYERECRSL TLREIHRFNN GLHSQNGYVT WNVDRLESAI RLGLNKVCEE GIRIDSIGID TWGVDFVLLD QQGQRVGLPV AYRDSRTNGL MAQAQQQLGK RDIYQRSGIQ FLPFNTIYQL RALTEQQPEL IPHIAHALLI PDYFSYRLTG KMNWEYTNAT TTQLVNINSD DWDESLLAWS GANKAWFGRP THPGNVIGHW ICPQGNEIPV VAVASHDTAS AVIASPLNGS RAAYLSSGTW SLMGFESQTP FTNDTALAAN ITNEGGAAGR YRVLKNIMGL WLLQRVLQER QINDLPALIA ATQALPACRF IINPNDDRFI NPDEMCSEIQ AACRETAQPI PESDAELARC IFDSLALLYA DVLHELAQLR GEDFSQLHIV GGGCQNALLN QLCADACGIR VIAGPVEAST LGNIGIQLMT LDELNNVDDF RQVVSTTANL TTFTPNPDSE IAHYVAQIHS TRQTKELCA
Source	Baculovirus
Target Names	rhaB
Protein Names	Recommended name: Rhamnulokinase EC= 2.7.1.5 Alternative name(s): Rhamnulose kinase
Expression Region	1-489
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.