



Recombinant Escherichia coli O6:K15:H31 3-dehydroquinate synthase (aroB)

Product Code	CSB-EP609781EGY-B
Abbreviation	aroB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TC78
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MERIVVTLGE RSYPIIASG LFNEPASFLP LKSGEQVMLV TNETLAPLYL DKVRGVLEQA GVNVDVILP DGEQYKSLAV LDTVFTALLQ KPHGRDITLV ALGGGVVDL TGFAAASYQR GVRFIQVPTT LLSQVDSSVG GKTAVNHPLG KNMIGAFYQP ASVVVDLDC LKTLPPRELAS GLAEVIKYGI ILDGAFNWL EENLDALLRL DGPAMAYCIR RCECELKPQVV AADERETGLR ALLNLGHTFG HAIEAEMGYG NWLHGEAVAA GMVMAARTSE RLGQFSSAET QRIITLLTRA GLPVNGPREM SAQAYLPHML RDKKVLAGE MRLILPLAIGK SEVRSGVSHE LVLNAIADCQ SA
Source	E.coli
Target Names	aroB
Protein Names	Recommended name: 3-dehydroquinate synthase EC= 4.2.3.4
Expression Region	1-362
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.