



Recombinant Escherichia coli O6:K15:H31 UPF0597 protein yhaM (yhaM)

Product Code	CSB-YP609791EGY
Abbreviation	yhaM
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TCZ4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MFDSTLNPLW QRYILAVQEE VKPALGCTEP ISLALAAAVA AAEELEGPVER VEAWVSPNLM KNGLGVTVP G TGMVGLPIAA ALGALGGNAN AGLEV LKDAT AQAI SDAKAL LAAGKVS VKI QEP CDEILFS RAKVWNGEKW ACVTIVGGHT NIVHIETHNG VVFTQQACVT EGEQESPLTV LSRTTLAEIL K FVNEVPFAA IRFILD SAKL NCALSQEGLS GNWGLHIGAT LEKQCARGLL AKDLSSSIVI RTSAASDARM GGATLPAMSN SGSGNQGITA TMPVVVVAEH FGADDERLAR ALMLSHLSAI YIHNQLPRLS ALCAATTAAM GAAAGMAWL V DGRYETISMA ISSMIGDVSG MICDGASNSC AMKVSTSASA AWKAVLMALD DTA VTGN EGI VAHDVEQSIA NLCALASHSM QQTDRQIIEI MASKAR
Source	Yeast
Target Names	yhaM
Protein Names	Recommended name: UPF0597 protein yhaM
Expression Region	1-436
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.