



Recombinant Escherichia coli O6:K15:H31 tRNA-modifying protein ygfZ (ygfZ)

Product Code	CSB-EP609797EGY
Abbreviation	ygfZ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TDV3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	AFTPFPPRQ PTASARLPLT LMTLDDWALA TITGADSEKY MQGQVTADVS QMTEDQHLLA AHCDKAGKMW SNLRLFRDGD GFAWIERRSV REPQLTELKK YAVFSKVTIA PDDERVLLGV AGFQARAALA NLFSELP SRE KQVVKEGATT LLWFEHPAER FLIVTDEATA NMLTDKLRGE AELNNSQQWL ALNIEAGFPV IDAANSQQFI PQATNLQALG GISFKKGCYT GQEMVARAKF RGANKRALWL LTGSASRLPE AGEDLELKMG ENWRRGTGTVL AAVKLEDGQV VVQVVMNNDM EPDSIFRVRD DANTLRIEPL PYSLEE
Source	E.coli
Target Names	ygfZ
Protein Names	Recommended name: tRNA-modifying protein ygfZ
Expression Region	2-326
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.