



# Recombinant Escherichia coli O6:K15:H31 Sulfate adenylyltransferase subunit 1 (cysN)

<b>Product Code</b>	CSB-EP609799EGY-B
<b>Abbreviation</b>	cysN
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TEA7
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MNTALAQQIA NEGGVEAWMI AQQHKSLRLF LTCGSVDDGK STLIGRLLHD TRQIYEDQLS SLHNSKRHG TQGEKLDLAL LVDGLQAERE QGITIDVAYR YFSTEKRKFI IADTPGHEQY TRNMATGAST CELAILLIDA RKGVLQDQTRR HSFISTLLGI KHLVVAINKM DLVDYSEETF TRIREDYLTFF AGQLPGNLDI RFVPLSALEG DNVASQSESM PWYSGPTLLE VLETVEIQRV VDAQPMRFPV QYVNRPNLDF RGYAGTLASG RVEVGQRVKV LPSGVESNVA RIVTFDGDRE EAFAGEAITL VLTDEIDISR GDLLLADEA LPAVQSASVD VVWMAEQPLS PGQSYDIKIA GKKTRARVDG IRYQVDINNL TQREVENLPL NGIGLVDLTF DEPLVLDRYQ QNPVTGGLIF IDRLSNVTVG AGMVHEPVSQ ATAAPSEFSA FELELNALVR RHFPHWGARD LLGDK
<b>Source</b>	E.coli
<b>Target Names</b>	cysN
<b>Protein Names</b>	Recommended name: Sulfate adenylyltransferase subunit 1 EC= 2.7.7.4 Alternative name(s): ATP-sulfurylase large subunit Sulfate adenylate transferase Short name= SAT
<b>Expression Region</b>	1-475
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.



Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.