



Recombinant Escherichia coli O6:K15:H31 Holo-[acyl-carrier-protein] synthase (acpS)

Product Code	CSB-BP609804EGY
Abbreviation	acpS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TES5
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MAILGLGTDI VEIARIEAVI ARSGERLARR VLSDNEWEIW KTHHQPVRFL AKRFAVKEAA AKAFGTGIRN GLAFNQFEVF NDELGKPRLR LWGEALKLAE KLGVVNMHVT LADERHYACA TVIIES
Source	Baculovirus
Target Names	acpS
Protein Names	Recommended name: Holo-[acyl-carrier-protein] synthase Short name= Holo-ACP synthase EC= 2.7.8.7 Alternative name(s): 4'-phosphopantetheinyl transferase AcpS
Expression Region	1-126
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.