



Recombinant Escherichia coli O6:K15:H31 UPF0304 protein yfbU (yfbU)

Product Code	CSB-MP609808EGY
Abbreviation	yfbU
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TFF2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MEMTNAQRLLSNQYKMMTM LDPANAERYR RLQTIERGY GLQMRELDRE FGELKEETCR TIIDIMEMYH ALHVSWSNLQ DQQSIDERRV TFLGFDAATE ARYLGYVRFM VNVEGRYTHF DAGTHGFNAQ TPMWEKYQRM LNVWHACPRQ YHLSANEINQ IINA
Source	Mammalian cell
Target Names	yfbU
Protein Names	Recommended name: UPF0304 protein yfbU
Expression Region	1-164
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.