



Recombinant Escherichia coli O6:K15:H31 Asparagine--tRNA ligase (asnS)

Product Code	CSB-BP609839EGY
Abbreviation	asnS
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TJC3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MSVVPVADVL QGRVAVDSEV TVRGWVRTRR DSKAGISFLA VYDGSCFDPV QAVINNSLPN YNEDVLRLLT GCSVIVTGKV VASPGQGQQF EIQASKVEVA GWVEDPDTYP MAAKRHSIEY LREVAHLRPR TNLIGAVARV RHTLAQALHR FFNEQGFFWV STPLITASDT EGAGEMFRVS TLDLENLPRN DQGKVDKDFDKD FFGKESFLTV SGQLNGETYA CALSKIYTFG PTFRAENSNT SRHLAEFWML EPEVAFANLN DIAGLAEAML KYVFKAVLEE RADDMKFFAE RVDKDAVSRL ERFIEADFAQ VDYTDAVTIL ENCGRKFENP VYWGVDLSSE HERYLAEHF KAPVVKNYP KDIKAFYMRL NEDGKTVAAM DVLAPGIGEI IGGSQREERL DVLDERMLEM GLNKEDYWWY RDLRRYGTVP HSGFGLGFER LIAYVTGVQN VRDVIPFPRT PRNASF
Source	Baculovirus
Target Names	asnS
Protein Names	Recommended name: Asparagine--tRNA ligase EC= 6.1.1.22 Alternative name(s): Asparaginyl-tRNA synthetase Short name= AsnRS
Expression Region	1-466
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.