



Recombinant Escherichia coli O6:K15:H31 Tetraacyldisaccharide 4'-kinase (lpxK)

Product Code	CSB-BP609840EGY
Abbreviation	lpxK
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TJD8
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MIEKIWSGES PLWRLLLPLS WLYGLVSGAI RLCYKCLKLR AWRAPIPVVV VGNLTAGGNG KTPVVVWLVE QLQQRGIRVG VVSRGYGGKA ESYPLLSAD TTTAQAGDEP VLIYQRTGAP VAVSPVRSDA VKAILAQHPD VQIIVTDDGL QHYRLARDVE IVAIDGVRRF GNGWWLPAGP MRERAGRLKS VDAVIVNGGV PRSGEIPMHL LPGQAVNLRG GTRCDVAQLE HVVAIAGIGH PPRFFATLKM CGVQPEKCVL LADHQLSLNHA DVSALVSAGQ TLMVTEKDAV KCRAFAEENW WYLPVDAQLS GDEPAKLLAQ LTSLASGH
Source	Baculovirus
Target Names	lpxK
Protein Names	Recommended name: Tetraacyldisaccharide 4'-kinase EC= 2.7.1.130 Alternative name(s): Lipid A 4'-kinase
Expression Region	1-328
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.