



Recombinant Escherichia coli O6:K15:H31 Apo-citrate lyase phosphoribosyl-dephospho-CoA transferase (citX)

Product Code	CSB-EP609845EGY-B
Abbreviation	citX
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TK58
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MHLLPELASH HAVSIPELLV SRDERQARQH AWLKRHPVPL VSFTVVAPGP IKDSEVTRRI FNHGVTALRA LAAKQGWIQ EQAALVSASG PEGMLSIAAP ARDLKLATIE LEHSHPLGRL WDIVLTPEG DILSRRDYSL PPRRCLLCEQ SAAVCARGKT HQLTDLLNRM EALLNDVDAC NVN
Source	E.coli
Target Names	citX
Protein Names	Recommended name: Apo-citrate lyase phosphoribosyl-dephospho-CoA transferase EC= 2.7.7.61 Alternative name(s): Apo-ACP nucleodityltransferase Holo-ACP synthase Holo-citrate lyase synthase
Expression Region	1-183
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.