



Recombinant Escherichia coli O6:K15:H31 Apo-citrate lyase phosphoribosyl-dephospho-CoA transferase (citX)

| | |
|--------------------------|---|
| Product Code | CSB-EP609845EGY |
| Abbreviation | citX |
| Storage | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |
| Uniprot No. | Q0TK58 |
| Product Type | Recombinant Protein |
| Immunogen Species | Escherichia coli O6:K15:H31 (strain 536 / UPEC) |
| Purity | >85% (SDS-PAGE) |
| Sequence | MHLLPELASH HAVSIPELLV SRDERQARQH AWLKRHPVPL VSFTVVAPGP IKDSEVTRRI FNHGVTALRA LAAKQGWIQ EQAALVSASG PEGMLSIAAP ARDLKLATIE LEHSHPLGRL WDIVLTPEG DILSRRDYSL PPRRCLLCEQ SAAVCARGKT HQLTDLLNRM EALLNDVDAC NVN |
| Source | E.coli |
| Target Names | citX |
| Protein Names | Recommended name: Apo-citrate lyase phosphoribosyl-dephospho-CoA transferase EC= 2.7.7.61 Alternative name(s): Apo-ACP nucleodityltransferase Holo-ACP synthase Holo-citrate lyase synthase |
| Expression Region | 1-183 |
| Notes | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week. |
| Tag Info | Tag type will be determined during the manufacturing process. |
| Protein Length | full length protein |
| Reconstitution | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference. |
| Shelf Life | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C. |