



# Recombinant Escherichia coli O6:K15:H31 Allantoinase (allB)

<b>Product Code</b>	CSB-YP609846EGY
<b>Abbreviation</b>	allB
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TKD0
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MSFDLIKNG TVILENEARV VDIKVDGKI AAIGQDLGDA KDVMASGLV VSPGMVDAHT HISEPGRSHW EGYETGTRAA AKGGITTMIE MPLNQLPATV DRASIELKFD AAKGKLTIDA AQLGGLVSYN IDRLHELDEV GVVGFKCFVA TCGDRGIDND FRDVNDWQFF KGAQKLGELG QPVLVHCENA LICDALGEEA KREGRVTAHD YVASRPVFTE VEAIRRVLYL AKVAGCRLHV CHVSSPEGVE EVTRARQEGQ DVTCESEPHY FVLDTDQFEE IGTAKCSPP IRDLNQKGM WEKLFNGEID CLVSDHSPCP PEMKAGNIMK AWGGIAGLQS CMDVMFDEAV QKRGMSLPMF GKLMTAAND IFGLQQKGR IAPGKDAFVF IQPNSSYVLT NDDLEYRHKV SPYVVRTIGA RITKTILRGD VIYDIEQGFP VAPKGQFILK HQQ
<b>Source</b>	Yeast
<b>Target Names</b>	allB
<b>Protein Names</b>	Recommended name: Allantoinase EC= 3.5.2.5 Alternative name(s): Allantoin-utilizing enzyme
<b>Expression Region</b>	1-453
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.