



Recombinant Escherichia coli O6:K15:H31 ATP-dependent Clp protease proteolytic subunit (clpP)

Product Code	CSB-MP609847EGY
Abbreviation	clpP
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TKK4
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	ALVPMV IEQTSRGERS FDIYSRLLKE RVIFLTGQVE DHMANLIVAQ MLFLEAENPE KDIYLYINSP GGVITAGMSI YDTMQFIKPD VSTICMGQAA SMGAFLLTAG AKGKRFLPN SRVMIHQPLG GYQQQATDIE IHAREILKVK GRMNELMALH TGQSLEQUIER DTERDRFLSA PEAVEYGLVD SILTHRN
Source	Mammalian cell
Target Names	clpP
Protein Names	Recommended name: ATP-dependent Clp protease proteolytic subunit EC=3.4.21.92 Alternative name(s): Endopeptidase Clp
Expression Region	15-207
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.