



Recombinant Escherichia coli O6:K15:H31 Thiamine import ATP-binding protein ThiQ (thiQ)

Product Code	CSB-BP609857EGY
Abbreviation	thiQ
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TLS2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	≥85% (SDS-PAGE)
Sequence	MLKLTDTWL YHHLPMRFSL TVERGEQVAI LGPSGAGKST LLNLIAGFLT PASGLLTIDD VDHTTTTPSR RPVSMLFQEN NLFSLTVAQ NIGLGLNPGL KLNAAQKKM HAIHQMGID NLMARLPGEL SGGQRQRVAL ARCLVREQPI LLLDEPFSAL DPALRQEMLT LVSTSCQQQK MTLMLVSHSV EDAARIATRS VVVADGRIAW QGKTDELLSG KASASALLGI KG
Source	Baculovirus
Target Names	thiQ
Protein Names	Recommended name: Thiamine import ATP-binding protein ThiQ EC= 3.6.3.-
Expression Region	1-232
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.