



# Recombinant Escherichia coli O6:K15:H31 Non-specific ribonucleoside hydrolase rihC (rihC)

<b>Product Code</b>	CSB-EP609860EGY-B
<b>Abbreviation</b>	rihC
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TLW1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MRLPIFLDTD PGIDDAVAIA AAIFTPELDL QLMTTVAGNV SVEKTRNAL QLLHFVNVDI PLAQGAHVPL VRAPRDAASV HGESGMAGYD FVEHNRQPLG IPAFLAIRDA LMRAPVPTL VAIGPLTNTA LLLSQCECK PYIRRLVIMG GSAGRGNCTP NAEFNIAADP EAAACVFRSG IEIVMCGLDV TNQAILTPDY LATLPELNRT GKMLHALFSH YRSGSMQSG LRMHDLCAIAW LVRPELFTLK PCFVAVETQG EFTSGTTVVD IDGCLGKPAN VQVALDLNVK GFQQWVAEVL ALVP
<b>Source</b>	E.coli
<b>Target Names</b>	rihC
<b>Protein Names</b>	Recommended name: Non-specific ribonucleoside hydrolase rihC EC= 3.2.-.- Alternative name(s): Purine/pyrimidine ribonucleoside hydrolase
<b>Expression Region</b>	1-304
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.