



Recombinant Human Protein ripply1 (RIPPLY1)

Product Code	CSB-YP609990HU
Abbreviation	RIPPLY1
Storage	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>
Uniprot No.	Q0D2K3
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	≥85% (SDS-PAGE)
Sequence	MDSAACAAAA TPVPALALAL APDLAQAPLA LPGLLSPSCL LSSGQEVNGS ERGTCLWRPW LSSTNDSRQ MRKLVDLAAG GATAAEVTKA ESKFHHPVRL FWPKSRSFDY LYSAGEILLQ NFPVQATINL YEDSDSEEEE EDEEQEDEEE K
Source	Yeast
Target Names	RIPPLY1
Protein Names	Recommended name: Protein ripply1
Expression Region	1-151
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	<p>The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.</p> <p>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.</p>