



Recombinant *Oryza sativa* subsp. japonica

Probable protein phosphatase 2C 62

(Os07g0507000, LOC_Os07g32380)

Product Code	CSB-MP609993OFG
Abbreviation	Os07g0507000, LOC_Os07g32380
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0D673
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MAGKEIYHKM KDKVKDAFSS SGPETGKGKT KLSGKRVKHG YHLVKGKSNH PMEDYLVAEY RQEGEHDGL FAIFDGHLGH TVPDFLRSHL FDNILKQPEF LSNPQAAIRN AYQLTDAKIL ESAAELGRGG STAVTAILIS SENSVNLVVA NVGDSRAVIS KSGVAKQLSV DHEPNKERHS IEKKGGFVSN LPGDVPRVDG QLAVARAFGD RSLKKHLSSE PDVVEEPIDE NTDFLILASD GLWKVMSNQE AVDEIKDFKD AQAAAKHLTE QAVNRKSKDD ISCIVVKFLC
Source	Mammalian cell
Target Names	Os07g0507000
Protein Names	Recommended name: Probable protein phosphatase 2C 62 Short name= OsPP2C62 EC= 3.1.3.16
Expression Region	1-290
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.