



Recombinant *Oryza sativa* subsp. japonica Zinc finger CCCH domain-containing protein 42 (Os06g0519400, LOC_Os06g32720)

Product Code	CSB-EP609996OFG-B
Abbreviation	Os06g0519400, LOC_Os06g32720
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0DBW8
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MMQILCECCN SDHRSSSTPM ATTTSSASD PAAISPTPSQ QHASSTVKTL DDRRPAGTSS SAGETEPKAA VEPQEYPRRP GVPDCSYIVE FGSCKFGMRC LYNHPAKHAG GCDKLEHPQR PGEHDCLHYL RFGRCKYGMN CRFNHPPDRL PQQQVYFPWK ACHCHHSEGK SEAEHVKLN LGLPLRPGTG LCSYIMNRGI CKFGSNCKFH HPNSGSGHEK WDGSLQTNQI SSGVNIYSVL DHGELNEQPV PSKDDFQVSF VQNIVGFNFY IWCDPELAP
Source	E.coli
Target Names	Os06g0519400
Protein Names	Recommended name: Zinc finger CCCH domain-containing protein 42 Short name= OsC3H42
Expression Region	1-279
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.