



# Recombinant *Oryza sativa* subsp. *japonica* Monothiol glutaredoxin-S11 (GRXS11)

|                          |                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
|--------------------------|---------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Code</b>      | CSB-EP610262OFG-B                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Abbreviation</b>      | GRXS11                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Storage</b>           | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.                                                                                                                                                                                                                                                                                                    |
| <b>Uniprot No.</b>       | Q0IWL9                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Product Type</b>      | Recombinant Protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Immunogen Species</b> | <i>Oryza sativa</i> subsp. <i>japonica</i> (Rice)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                     |
| <b>Purity</b>            | >85% (SDS-PAGE)                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                       |
| <b>Sequence</b>          | MAAVREVGSK AELEAAAGGA RAAAVHFWAA WCEASKQMDE<br>VFAHLAVDFS HAVFLRVEAE EQPEISEAYG VTAVPYFVFL KEGKTVDTLE<br>GANPASLANK VAKLAGPASV AESAVPASLG VAAGPAVLEK VQEMAQQNGA<br>SATSSAEDAL NKRLEQLVNS HPVFLFMKGT PEQPRCGFSR KVVDVLKQEG<br>VEFGSFDILT DNDVREGMCK FSNWPTFPQL YCKGELLGGC DIVIAMHESG<br>ELKDVFKEHN IPLQPQGSKN EEAVKAKPDT EKSGAVSEPA LLTAAQKERL<br>ESLVNFSTVM AFIKGTPEEP KCGFSGKLVH ILKQEKIPFS SFDILTDDEV<br>RQGLKLLSNW PSYPQLYING ELVGGSDIVM EMHKSSELK VLSEKGVAK<br>ESLEDRLKAL ISSAPVMLFM KGTPDAPRCG FSSKVVNALK QAGVSFGAFD<br>ILSDEEVRQG LKTYSNWPTF PQLYYKSELI GGCDIVLELE KSGELKSTLS E |
| <b>Source</b>            | E.coli                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Target Names</b>      | GRXS11                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                |
| <b>Protein Names</b>     | Recommended name: Monothiol glutaredoxin-S11                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                          |
| <b>Expression Region</b> | 1-491                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                 |
| <b>Notes</b>             | Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Tag Info</b>          | Tag type will be determined during the manufacturing process.                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                         |
| <b>Protein Length</b>    | full length protein                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                                   |
| <b>Reconstitution</b>    | We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.                                                                                                                                                                           |
| <b>Shelf Life</b>        | The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.<br>Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life                                                                                                                                                                                                                                                                                                                                                     |



of lyophilized form is 12 months at -20°C/-80°C.