



Recombinant *Oryza sativa* subsp. japonica Coatomer subunit delta-4 (Os08g0368000, LOC_Os08g28080)

Product Code	CSB-BP610264OFG
Abbreviation	Os08g0368000, LOC_Os08g28080
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0J649
Product Type	Recombinant Protein
Immunogen Species	<i>Oryza sativa</i> subsp. japonica (Rice)
Purity	>85% (SDS-PAGE)
Sequence	MTFLYITCYC VAFNRKLLVF MQMEKGKHNK GGYSSISGPQ VIEKSFNDMS ISCTGFGSGS GLGGLNTDTD TFTSRPKGRT SGGTTGAGKG IGMKLGNTKK TNQFLESLKA EGEVIMEDFQ PCSLQSRSSP LPPSDPVTVA VEEKLNVAVK RDGGVNNFDV QGTLALLVLN DADGLILLQI ESQDIPGLSF KTHPNINKDL FNSQQILGAK DPIRPFPSGQ NETPLVKWRI QGMNESSLPL SVNCWPSILG NETYVNIEYE ASEMFDLHSV IISIPLPALR EAPRVRQIDG EWKYDSRNSV LEWSIILIDQ SNRSGSMEFV VPPADPSMFY PISVGFASAN TFSNVKVTGI RPLNEGSNPP KYSQRVRLVA DNYQVV
Source	Baculovirus
Target Names	Os08g0368000
Protein Names	Recommended name: Coatomer subunit delta-4 Alternative name(s): Delta-coat protein 4 Short name= Delta-COP 4
Expression Region	1-376
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.