



Recombinant Human Protein ABHD18 (ABHD18)

Product Code	CSB-BP610370HU
Abbreviation	C4orf29
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0P651
Product Type	Recombinant Protein
Immunogen Species	Homo sapiens (Human)
Purity	>85% (SDS-PAGE)
Sequence	RPEDLK RLFEFRKMIG NRERCQNLVS SDYPVHIDKI EEQSDCKILD GHFVSPMAHY VPDIMPIESV IARFQFIVPK EWNSKYRPVC IHLAGTGDHH YWRRRTL MAR PMIKEARMAS LLENPPYYGC RPKKDQVRSS LKNVSDLFVM GGALVLESAA LLHWLREGY GPLGMTGISM GGHMASLAVS NWPKPMPLIP CLSWSTASGV FTTTDSFKMG QEFVKHFTSS ADKLTNLNLV SRTLNLDISN QVVSQKPADC HNSSKTSVSA TSEGLLLQDT SKMKRFNQTL STNKSGYTSR NPQSYHLLSK EQSRNSLRKE SLIFMKGVM D ECTHVANFSV PVDPSLIIVV QAKEDAYIPR TGVRSLQEIW PGCEIRYLEG GHISAYLFKQ GLFR
Source	Baculovirus
Target Names	ABHD18
Protein Names	Recommended name: Uncharacterized protein C4orf29
Expression Region	25-414
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.