



Recombinant Escherichia coli O6:K15:H31 Glucose-1-phosphate adenylyltransferase (glgC)

Product Code	CSB-MP610637EGY
Abbreviation	glgC
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TC29
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MVSLEKNDHL MLARQLPLKS VALILAGGRG TRLKDLTNKR AKPAVHFGGK FRIIDFALSN CINSGIRRMG VITQYQSHTL VQHIQRGWSF FNEEMNEFVD LLPAQQRMKG ENWYRGTADA VTQNLDIIRR YKAEYVVILA GDHIYKQDYS RMLIDHVEKG ARCTVACMPV PIEEASAFGV MAVDENDKII EFVEKPANPP SMPNDPSKSL ASMGYVFDA DYLYELLEED DRDENSSHDF GKDLIPKITE AGLAYAHPFP LSCVQSDPDA EPYWRDVGTL EAYWKANLDL ASVVPEDMY DRNWPIRTYN ESLPPAKFVQ DRSGSHGMTL NSLVSGGCVI SGSVVVQSVL FSRVRVNSFC NIDSAVLLPE VVWGRSCLLR RCVIDRACVI PEGMVIGENA EEDARRFYRS EEGIVLVTRE MLRKLGHKQE R
Source	Mammalian cell
Target Names	glgC
Protein Names	Recommended name: Glucose-1-phosphate adenylyltransferase EC= 2.7.7.27 Alternative name(s): ADP-glucose pyrophosphorylase Short name= ADPGlc PPase ADP-glucose synthase
Expression Region	1-431
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.