



Recombinant Escherichia coli O6:K15:H31 Ribosomal RNA small subunit methyltransferase B (rsmB)

Product Code	CSB-YP610644EGY
Abbreviation	rsmB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TCH3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MKKQRNLRSM AAQAIEQVVE QGQSLSNILP PLQQKVSDKD KALLQELCFG VLRTLSQLDW LINKLMARPM TGKQRTVHYL IMVGLYQLLY TRIPPHAALA ETVEGAVAIK RPQLKGLING VLRQFQRQQD ELLAEFNASD ARYLHPSWLL KRLQKAYPEQ WQSIVEANNQ RPPMWLRVNR THHSRDSWLA LLDEAGMKGF PHADYPDAVQ LETPAPVHAL PGFEEGWVTV QDASAQGCMT WLAPQNGEHI LDLCAGPGGK TTHILEVAPE AQVLAVDIDE QRLSRVYDNL KRLGMKATVK QGDGRYPSQW CGEQQFDRIL LDAPCSATGV IRRHPDIKWL RRDRDIPELA QLQSEILDAI WSHLKSGGTL VYATCSVLPE ENSLQIKAFI QRTADAELCE TGTPEQPGKQ NLPGAEEGDG FFYAKLIKK
Source	Yeast
Target Names	rsmB
Protein Names	Recommended name: Ribosomal RNA small subunit methyltransferase B EC=2.1.1.176 Alternative name(s): 16S rRNA m5C967 methyltransferase rRNA (cytosine-C(5)-)-methyltransferase RsmB
Expression Region	1-429
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.