



Recombinant Escherichia coli O6:K15:H31 UDP-N-acetylglucosamine 1-carboxyvinyltransferase (murA)

Product Code	CSB-BP610647EGY
Abbreviation	murA
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TCS1
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MDKFRVQGPT KLQGEVTISG AKNAALPILF AALLAEEPVE IQNVPKLDV DTSMKLLSQL GAKVERNCSV HIDARDVNVF CAPYDLVKTM RASIWALGPL VARFGQQQVS LPGGCTIGAR PVDLHISGLE QLGATIKLEE GYVKASVDGR LKGAHIVMDK VSVGATVTIM CAATLAEGTT IIENAAREPE IVDTANFLIT LGAKISGQGT DRIVIEGVER LGGGVYRVLP DRIETGFTLV AAISRKII CRNAQPDTLD AVLAKLRDAG ADIEVGEDWI SLDMHGKRPK AVNVRTAPHP AFPTDMAQAF TLLNLVAEGT GFITETVFEN RFMHVPELSR MGAHAEIESN TVICHGVEKL SGAQVMATDL RASASLVLAG CIAEGTTVVD RIYHIDRGYE RIEDKLRALG ANIERVKGE
Source	Baculovirus
Target Names	murA
Protein Names	Recommended name: UDP-N-acetylglucosamine 1-carboxyvinyltransferase EC= 2.5.1.7 Alternative name(s): Enolpyruvate transferase UDP-N- acetylglucosamine enolpyruvyl transferase Short name= EPT
Expression Region	1-419
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.