



# Recombinant Escherichia coli O6:K15:H31 UDP-N-acetylglucosamine 1-carboxyvinyltransferase (murA)

<b>Product Code</b>	CSB-MP610647EGY
<b>Abbreviation</b>	murA
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TCS1
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	MDKFRVQGPT KLQGEVTISG AKNAALPILF AALLAEEPVE IQNVPKLDV DTSMKLLSQL GAKVERNCSV HIDARDVNVF CAPYDLVKTM RASIWALGPL VARFGQQQVS LPGGCTIGAR PVDLHISGLE QLGATIKLEE GYVKASVDGR LKGAHIVMDK VSVGATVTIM CAATLAEGTT IIENAAREPE IVDTANFLIT LGAKISGQGT DRIVIEGVER LGGGVYRVLP DRIETGTLV AAAISRGKII CRNAQPDTLD AVLAKLRDAG ADIEVGEDWI SLDMHGKRPK AVNVRTAPHP AFPTDMAQAF TLLNLVAEGT GFITETVFEN RFMHVPELSR MGAHAEIESN TVICHGVEKL SGAQVMATDL RASASLVLAG CIAEGTTVVD RIYHIDRGYE RIEDKLRALG ANIERVKGE
<b>Source</b>	Mammalian cell
<b>Target Names</b>	murA
<b>Protein Names</b>	Recommended name: UDP-N-acetylglucosamine 1-carboxyvinyltransferase EC= 2.5.1.7 Alternative name(s): Enolpyruvate transferase UDP-N- acetylglucosamine enolpyruvyl transferase Short name= EPT
<b>Expression Region</b>	1-419
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	full length protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.



## Shelf Life

The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself.

Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.