



Recombinant Escherichia coli O6:K15:H31 Isopentenyl-diphosphate Delta-isomerase (idi)

Product Code	CSB-EP610656EGY-B
Abbreviation	idi
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TDW3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MQTEHVILLN AQQVPTGTLE KYAAHTADTL LHAFSSWLF NAKGQLLVTR RALSKKAWPG MWTNSVCGHP QLGESNEEAV IRRCRYELGV EITPPESIYP DFRYRATDPN GIVENEVCPV FAARTTSALQ INDDEVMDYQ WCDLAAVLRG IDATPWAFSP WMVMQATNRE ARKRLSAFTQ LK
Source	E.coli
Target Names	idi
Protein Names	Recommended name: Isopentenyl-diphosphate Delta-isomerase Short name= IPP isomerase EC= 5.3.3.2 Alternative name(s): IPP:DMAPP isomerase Isopentenyl pyrophosphate isomerase
Expression Region	1-182
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.