



Recombinant Escherichia coli O6:K15:H31 Probable transcriptional regulatory protein YeeN (yeeN)

Product Code	CSB-EP610677EGY
Abbreviation	yeeN
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TGH3
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MGRKWANIVA KKTAKDGATS KİYAKFGVEI YAAAKQGEPD PELNTSLKFV IERAKQAQVP KHVIDKAIDK AKGGGDETFV QGRYEGFGPN GSMIIAETLT SNVNRTIANV RTIFNKKGGN IGAAGSVSYM FDNTGVIVFK GTDPDHIFEI LLEAEVDVRD VTEEEGNIVI YTEPTDLHKG IAALKAAGIT EFSTTELEMI AQSEVELSPE DLEIFEGLVD ALEDDDDVQK VYHNVANL
Source	E.coli
Target Names	yeeN
Protein Names	Recommended name: Probable transcriptional regulatory protein YeeN
Expression Region	1-238
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.