



Recombinant Escherichia coli O6:K15:H31 8-amino-7-oxononanoate synthase (bioF)

Product Code	CSB-BP610700EGY
Abbreviation	bioF
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TJS2
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MIWQEKIDAA LDARRVADAL RRRYPVAQGA GRWLVADDCQ YLNFSSNDYL GLSHHPQIIR AWQQGAEQFG VGSGGSGHVS GYSVAHQVLE EELAEWLGYS RALLFISGFA ANQAVIAAMM AKEDRIVADR LSHASLLEAA SLSPSPLRRF AHNDVTHLAR LLASPCPGQQ LVVTEGVFSM DGDSAPLEEI QQVTQQHGDGW LMVDDAHGTG VIGEQQRGSC WLQKVKPELL VVTFGKGFV SGAAVLCSNT VADYLLQFAR HLIYSTSMPP AQAQALRASL AVIRSDEGDA RREKLAALIT RFRAGVQDLP FTLADSWSAI QPLIVGDNRSR ALQLAEKLRQ QGCWVTAIRP PTVPAGTARL RLTLTAAHEM QDIDRLLEVL HGNG
Source	Baculovirus
Target Names	bioF
Protein Names	Recommended name: 8-amino-7-oxononanoate synthase Short name= AONS EC= 2.3.1.47 Alternative name(s): 7-keto-8-amino-pelargonic acid synthase Short name= 7-KAP synthase Short name= KAPA synthase 8-amino-7- ketopelargo
Expression Region	1-384
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.