



Recombinant Escherichia coli O6:K15:H31 Ribonuclease HII (rnhB)

Product Code	CSB-MP610712EGY
Abbreviation	rnhB
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0TLF0
Product Type	Recombinant Protein
Immunogen Species	Escherichia coli O6:K15:H31 (strain 536 / UPEC)
Purity	>85% (SDS-PAGE)
Sequence	MIEFVYPHTQ LVAGVDEVGR GPLVGAVVTA AVILDPARPI AGLNDSKKLS EKRRALCEE IKEKALSWSL GRAEPHEIDE LNILHATMLA MQRVAVAGLHI APEYVLIDGN RCPKLPMP SM AVVKGDSRVP EISAASILAK VTRDAEMAAL DIVFPQYGFA QHKGYP TAFH LEKLAEHGAT EHHRRSFGPV KRALGLAS
Source	Mammalian cell
Target Names	rnhB
Protein Names	Recommended name: Ribonuclease HII Short name= RNase HII EC= 3.1.26.4
Expression Region	1-198
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	full length protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.