



# Recombinant Clostridium perfringens Phospholipase C (plc)

<b>Product Code</b>	CSB-EP610762CAAN-B
<b>Abbreviation</b>	plc
<b>Storage</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
<b>Uniprot No.</b>	Q0TV31
<b>Product Type</b>	Recombinant Protein
<b>Immunogen Species</b>	Clostridium perfringens (strain ATCC 13124 / DSM 756 / JCM 1290 / NCIMB 6125 / NCTC 8237 / Type A)
<b>Purity</b>	>85% (SDS-PAGE)
<b>Sequence</b>	WD GKIDGTGTHA MIVTQGVSI L ENDSLKNEPE SVRKNLEILK ENMHELQLGS TYPDYDKNAY DLYQDHFWD P DTDNNFSKDN SWYLAYSIPD TGESQIRKFS ALARYEWQRG NYKQATFYLG EAMHYFGDID TPYHPANVTA VDSAGHVKFE TFAEERKEQY KINTAGCKTN EAFYTDILKN KDFNAWSKEY ARGFAKTGKS IYSHASMSH SWDDWDYAAK VTLANSQKGT AGYIYRFLHD VSEGNDPSVG KNVKELVAYI STSGEKDAGT DDYMYFGIKT KDGKTQEWEM DNPGNDFMTG SKDITYFKLK DENLKIDDIQ NMWIRKRKYT AFSDAYKPEN IKIIANGKVV VDKDINEWIS GNSTYNIK
<b>Source</b>	E.coli
<b>Target Names</b>	plc
<b>Protein Names</b>	Recommended name: Phospholipase C Short name= PLC EC= 3.1.4.3 Alternative name(s): Alpha-toxin Hemolysin Lecithinase Phosphatidylcholine cholinephosphohydrolase
<b>Expression Region</b>	29-398
<b>Notes</b>	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
<b>Tag Info</b>	Tag type will be determined during the manufacturing process.
<b>Protein Length</b>	Full Length of Mature Protein
<b>Reconstitution</b>	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
<b>Shelf Life</b>	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life



of lyophilized form is 12 months at -20°C/-80°C.