



Recombinant Arabidopsis thaliana Protein FAF-like, chloroplastic (At5g22090)

Product Code	CSB-EP610794DOA-B
Abbreviation	At5g22090
Storage	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.
Uniprot No.	Q0V865
Product Type	Recombinant Protein
Immunogen Species	Arabidopsis thaliana (Mouse-ear cress)
Purity	>85% (SDS-PAGE)
Sequence	ADLSSK TWVSQNGFSP MKRISSEKL RPDEEEAEEE SRSGVDIWAQ IQQDKNDKKK EEEIEPGQSD VWSSILSEKK KTESSKDTVP PPYVHPLMKR ASSLSEKSLE ICTESLGSET GCDGFSSHAS SETGDAEIEI HDETNLVVNV TETKVEEITE TEIVVEQESS IIVPNHIEL PPGSFPPPIR SLSSQSGSSL HMKTRRDNGR LVLEAVSMPS HNNFSAKRQD GRLLTFAEI SNEPNYDKED EIDSEVQWFD EEEEEEEEEEE DEEEEAPDEF AYKPNGLLYK MAQKPIPIV HRLAHKPIGV PKRNSRWPMA DEFDTKSDLS TPVVHSLPPR PRVAQLARST KPPSTVDDTV GAACFNTCDY SWKSTNTESF GPNTKTQFQA QNFVNKSMGD GWINGCKDRR RLLSVEPFC IAT
Source	E.coli
Target Names	At5g22090
Protein Names	Recommended name: Protein FAF-like, chloroplastic
Expression Region	45-463
Notes	Repeated freezing and thawing is not recommended. Store working aliquots at 4°C for up to one week.
Tag Info	Tag type will be determined during the manufacturing process.
Protein Length	Full Length of Mature Protein
Reconstitution	We recommend that this vial be briefly centrifuged prior to opening to bring the contents to the bottom. Please reconstitute protein in deionized sterile water to a concentration of 0.1-1.0 mg/mL. We recommend to add 5-50% of glycerol (final concentration) and aliquot for long-term storage at -20°C/-80°C. Our default final concentration of glycerol is 50%. Customers could use it as reference.
Shelf Life	The shelf life is related to many factors, storage state, buffer ingredients, storage temperature and the stability of the protein itself. Generally, the shelf life of liquid form is 6 months at -20°C/-80°C. The shelf life of lyophilized form is 12 months at -20°C/-80°C.